

POSITIONING SUPPORTS



Superior Clinical Seating.

ADJUSTABLE STRETCH, CENTER-OPENING



	DYNAMIC	STRUCTURED	A	B
EI	ATS8807	ATS9407	5"	9.75"
XS	ATS8800	ATS9400	5.75"	11"
S	ATS8801	ATS9401	6.5"	12.25"
M	ATS8802	ATS9402	9.75"	14.25"
L	ATS8803	ATS9403	9.75"	15.5"
XL	ATS8804	ATS9404	10"	17.5"

The Center Opening allows for increased flexibility in frequent transfers and the adjustable center strap allows for a more precise fit.

TRUNK STRAP

The Trunk Strap style allows for increased range of motion when shoulder rotation and reaching are required.



	STANDARD	BETWEEN THORACICS	A	B
XS	ATS4050	ATS4050T	8"	3"
S	ATS4000	ATS4000T	10"	3"
M	ATS4100	ATS4100T	12"	3"
L	ATS4200	ATS4200T	14"	4"
XL	ATS4300	ATS4300T	16"	4"
Custom	ATS4900	ATS4900T	<i>Specify</i>	

BACKPACK STYLE

Allows control of forward flexion and allows more intimate shoulder retraction—can be used aggressively for upper torso control.



	DYNAMIC	STRUCTURED	A
S	ATS9601	ATS9701	12"
L	ATS9603	ATS9703	18"

Y-STYLE



	DYNAMIC	STRUCTURED	A	B
EI	ATS9167	ATS9367	15"	7"
S	ATS9161	ATS9361	19"	8"
M	ATS9162	ATS9362	22"	9"
L	ATS9163	ATS9363	25"	11"
Custom	ATS9169	ATS9369	<i>Specify</i>	

The Y-Style allows for increased range of motion for one arm while still providing optimal support.

Hand and Arm Supports

Hand & Arm Supports help maintain a comfortable and optimal position of the hands and arms for improved upper body posture and head control.

ARM TROUGHS



Fixed Hardware



Swing-Away Hardware

A	HARDWARE	
	FIXED	SWING-AWAY
14"	JARM14F	JARM14SA
18"	JARM18F	JARM18SA
22"	JARM22F	JARM22SA

PADDED WRIST STRAPS



	STRAP	WRIST CIRCUMFERENCE
S	ATS7510	5.5" to 7"
M	ATS7520	7.5" to 9"
L	ATS7530	9.5" to 11"
Custom	ATS7590	<i>Specify</i>

GEL ARMREST PAD COVERS



	GEL PAD COVER
DESK	131066
FULL	131067
PWR RECLINE (R)	131071
PWR RECLINE (L)	131072

LOWER BODY POSITIONING



Shoe Holders

Shoe Holders help maintain positioning of the feet and ankles, broadening the seating "footprint", and improving overall posture.

FOOT AND ANKLE POSITIONING (EACH)



	FOOT & ANKLE	ANKLE ONLY	A	B
XS	—	TST6501R TST6501L	6"	4.5"
S	TST6510R TST6510L	TST6511R TST6511L	7.75"	5"
M	TST6520R TST6520L	TST6521R TST6521L	7.75"	5"
L	TST6540R TST6540L	TST6541R TST6541L	9.75"	6"
XL	TST6560R TST6560L	TST6561R TST6561L	10.75"	6.5"

SHOE HOLDER (EACH)



	PADDED	NON-PADDED	A	B
S	TST6001	TST6000	5.75"	3"
M	TST6101	TST6100	7.75"	3.25"
L	TST6201	TST6200	9.75"	4"
XL	TST6301	TST6300	12.75"	6"

FLEXI-FEET (PAIR)



	FLEXI-FEET	A	B
XS	32600XS-RT	4.5"	2.75"
S	32600S-RT	6.25"	3"
M	32600M-RT	8.25"	3.5"
L	32600L-RT	9"	4.5"

Foot Boxes

Foot Boxes reduce the risk of skin breakdown and/or client injury that may result from periods of high tone or continual contact with the hangers and footplates.

SINGLE FOOT BOX



CALCULATING FOOT BOX WIDTH	
Flared Frame	Foot Box Width = Seat Width + 2"
Offset Hanger	Foot Box Width = Seat Width + 2"
0" Inset	Foot Box Width = Seat Width
1" Inset	Foot Box Width = Seat Width - 2"
2" Inset	Foot Box Width = Seat Width - 4"

	HEIGHT	A	B
JFOOTSS	Standard	8"	7.5"
JFOOTST	Tall	12"	7.5"

FULL FOOT BOX W/DIVIDER



	WIDTH	HEIGHT	A	B
JFOOT16SD	16"	Standard	8"	13.5"
JFOOT18SD	18"	Standard	8"	15.5"
JFOOT20SD	20"	Standard	8"	17.5"
JFOOT22SD	22"	Standard	8"	19.5"
JFOOT16TD	16"	Tall	12"	13.5"
JFOOT18TD	18"	Tall	12"	15.5"
JFOOT20TD	20"	Tall	12"	17.5"
JFOOT22TD	22"	Tall	12"	19.5"

FULL FOOT BOX



	WIDTH	HEIGHT	A	B
JFOOT12S	12"	Standard	8"	9.5"
JFOOT14S	14"	Standard	8"	11.5"
JFOOT16S	16"	Standard	8"	13.5"
JFOOT18S	18"	Standard	8"	15.5"
JFOOT20S	20"	Standard	8"	17.5"
JFOOT22S	22"	Standard	8"	19.5"
JFOOT16T	16"	Tall	12"	13.5"
JFOOT18T	18"	Tall	12"	15.5"
JFOOT20T	20"	Tall	12"	17.5"
JFOOT22T	22"	Tall	12"	19.5"